ELANEE®

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## ELANEE Beckenboden-Trainingshilfen

### Phase I aktiv

<b>DE</b>	Beckenboden-Trainingshilfe Phase I aktiv – Gebrauchsanweisung.....	2
<b>EN</b>	Pelvic Floor Training Aids Phase I active – Instructions for Use.....	9
<b>FR</b>	Accessoires de rééducation du plancher pelvien Niveau I actif – Notice d'utilisation.....	16
<b>IT</b>	Ausili per l'allenamento dei muscoli del pavimento pelvico fase I attiva – Istruzioni per l'uso.....	24
<b>NL</b>	Bekkenbodem-trainingshulpjes Fase I actief – Gebruiksaanwijzing.....	32
<b>ES</b>	Ayudas de entrenamiento del suelo pélvico - fase I activa – Instrucciones de uso.....	40
<b>RU</b>	Конусы для укрепления мышц тазового дна Фаза I активная – Инструкция по применен.....	48

**LOT**  siehe Verpackung / see packaging / voir l'emballage / vedere  
imballaggi / zie verpakking / Véase el embalaje / См. упаковку

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Artikelnummer

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---

## **EN Pelvic Floor Training Aids Phase I active – Instructions for use**

**Dear Customer,**

thank you for choosing our Pelvic Floor Training Aids Phase I active. The Instructions for Use should help you to get the best

use out of your product. Please read the instruction manual carefully – particularly the usage and safety instructions – before using this product for the first time.

We will be glad to answer any questions and to receive any suggestions you might have

*Your ELANEE-Team*

 Please follow the instructions and retain the instruction sheet and the packaging!

*Made in Germany*

### **1. Field of application**

These instructions for use apply to the following product:

**REF** 701-00/702-00 ELANEE Pelvic Floor Training Aids Phase I active

### **2. Product content**

ELANEE Pelvic Floor Training Aids Phase I active consist of a set of four tampon-shaped cones with a different colour and weight in each case:

<b>Colour</b>	yellow	blue	green	violet
<b>Weight</b>	20g	28g	52g	71g

In addition, Instructions for Use and an exercise booklet are enclosed.

### **3. Purpose and intended use**

ELANEE Phase I active pelvic floor training aids are a set of four equal-sized tampon-shaped cones with different weights to be used in the vagina for short periods of time.

The training is for the prevention, reduction (alleviation) or correction (treatment) of stress incontinence, caused by weak muscles in the pelvic floor.

The training is suitable for women

- during pregnancy and after giving birth to prevent pelvic floor weakness and to support regeneration and/or
- to prevent, reduce or correct bladder weakness due to stress incontinence, in particular during pre-menopause or the menopause.

The product is only intended for private home use. The ELANEE Phase I active pelvic floor training aids are not suitable for any purposes other than the listed intended ones.

#### 4. Areas of application

- **During an unproblematic pregnancy** the training with ELANEE Pelvic Floor Training Aids Phase I active can strengthen the pelvic floor and prevent potential subsequent problems.
- **After giving birth you can resume the training as soon as the cervix is closed (after about 6 week postpartum period).** If in doubt, ask your doctor or midwife.
- Training with ELANEE Pelvic Floor Training Aids Phase I active can bring about the desired success in cases of **ineffective non-invasive pelvic floor training** caused by not being able to feel the pelvic floor muscles.
- **The ELANEE Pelvic Floor Training Aids Phase I active are also suitable for the treatment of weaknesses of the pelvic floor resulting from a weak bladder (load/stress**

**incontinence).** However, bladder problems have various causes. Other causes of the existing urinary incontinence should be ruled out by your doctor and/or midwife in advance.

## **5. Functionality**

The pelvic floor is a flat, elastic muscle plate which is approximately palm-sized located between the pubic bone, coccyx and the buttock. It supports the internal organs in the abdominal cavity. The urethra, vagina and anus pass through the pelvic. Pregnancy and childbirth, hormonal changes/deficiency during the climacteric period or menopause, heavy physical exertion or being overweight can lead to over-stretching or weakening of the pelvic floor. The consequences can be urinary incontinence (load/ stress incontinence) or an incorrect position of the organs in the abdomen. However, the pelvic floor can be strengthened through exercise, just like any other muscle. A strong pelvic floor can also have a positive impact on sexual pleasure.

The ELANEE Phase I active pelvic floor training aids consist of a set of four different weights, which enable active pelvic floor training in accordance with the biofeedback principle.

Underlying the biofeedback principle is the fact that the cones will slide out of the vagina due to their own weight. The muscles of the pelvic floor react actively to this, by spontaneously and automatically contracting even tighter, in order to hold the weight. These repeated active contractions ensure muscle development – just like when muscle training with barbells.

Four different weights in the kit enable you to increase active exercise intensity. The training weight is increased if the muscles can hold a cone with ease. In this way a huge increase can

be achieved from the lightest to the heaviest weight. An active training with pelvic floor training aids makes sense for women of all ages. It is used during pregnancy and after birth to prevent weaknesses of the pelvic floor and facilitate regeneration. In addition, it is ideal for reducing, eliminating and preventing urinary incontinence (due to load/stress incontinence) in women, particularly during the climacteric period and menopause.

## 6. Usage

1. Thoroughly **wash** the ELANEE Pelvic Floor Training Aids prior to initial usage and ensure clean conditions (e.g. wash hands).
2. **Test** to see which cone you will use to start your workout: To do this, take the **yellow (lightest) cone** and insert it slowly, completely and deeply about 3-4 cm into the vagina as you would normally insert a tampon.
  - If you can hold the yellow cone for 10 minutes while standing and walking, you can use the next heaviest cone (blue).
  - If you can also hold the blue cone for 10 minutes, take the next weight, etc.
3. You complete the training with the cones by merely increasing the carrying time until the cone being used can be held **twice daily for over 10 minutes** while standing and walking, for several days in succession. Now consciously increase your feeling of your pelvic muscles by actively contracting and relaxing your muscles in alternation in order to let the cone slide out and draw in again.
4. **After a few weeks** of daily training with the ELANEE Pelvic Floor Training Aids Phase I active you will have learnt to

feel the pelvic floor muscles without the use of cones. However, occasional practice with the help of the cones is still recommended.

**Tip:** The use of medical lubricating gel can be helpful when training with the ELANEE Phase I active pelvic floor training aids, in particular if your vagina is dry.

## 7. Safety Instructions

-  Thoroughly clean the ELANEE Pelvic Floor Training Aids Phase I activeboth prior to initial usage and after every application and ensure clean conditions (e.g. wash hands).
-  Please check your training aids before each use. Only use cones which are in perfect condition. A training aid with a scratched or damaged surface could cause injuries.
-  In cases of vaginal infections, inflammations, injury or painful complaints in the genital area and during menstruation or a problematic pregnancy, you should avoid using the Pelvic Floor Training Aids.
-  When other complaints (such as a dropping of the pelvic organs) exist, you should have other causes than a weakness of the pelvic floor ruled out in advance by a doctor before you start with the training.
-  In cases of a known hypersensitivity and allergies to one of the named ingredients, you should not use the aids or only use them after consultation with a doctor.
-  If you are a virgin, you should not use the pelvic floor training aids.
-  The wearing of cervical caps, contraceptive rings, dia-

phragms or intrauterine devices during training are possible. However, you must check the correct positioning before and after each training session (if in doubt only use after consultation with a doctor).

-  You must not use the cones during X-rays, a CT scan (computer tomography) or an MRI scan (magnetic resonance imaging).
-  If the symptoms persist or worsen, please seek advice from a doctor.
-  All serious incidents occurring in connection with this product must be reported to the manufacturer and the responsible authorities in the user's country of residence.
-  Do not use steam sterilisation, dry heat sterilisation or boiling to clean the aids as this can lead to distortion of the material.

## **8. Safekeeping, service life, storage and transport**

There are no special instructions for keeping, storing and transporting the aids. With normal usage, the training aids should be replaced after 5 years.

## **9. Cleaning and disposal instructions**

Clean the training aids after each use with warm water and soap or disinfect the cone with a commercially available skin disinfectant.

Before being used again, any possible soap residue or disinfectants should be thoroughly rinsed off the training aids with water.

The ELANEE Pelvic Floor Training Aids Phase I active can be disposed of as standard domestic waste.

Please separate the individual components of the packaging into groups of identical materials and dispose of them in the appropriate recycling system.

When doing this please note the local and regional waste disposal regulations.

## 10. Material composition

The ELANEE Pelvic Floor Training Aids Phase I active consist of:

- an aluminium or steel core coated with polypropylene
- a polyamide yarn coated with polyurethane

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